

BIONOMADIX

*Wireless, Wearable
Physiology Monitoring*

C A T A L O G

For the Life Sciences



- FULL-BANDWIDTH, HIGH-QUALITY DATA FOR A WIDE RANGE OF SIGNALS
- COMFORTABLE FOR THE SUBJECT AND EMPOWERING FOR THE RESEARCHER
- USE FOR IN-LAB OR REAL LIFE STUDIES, INCLUDING ACTIVE OR LONG-TERM PROTOCOLS

**NEW SYSTEMS FOR
IN LAB, PORTABLE,
OR LOGGED DATA**

BIOPAC[®]
Systems, Inc.
Registered to ISO 9001:2015



EXPERIENCE THE FREEDOM OF DISCOVERY



Linga © gert. weigel!

Advanced Wireless, Wearable Solutions for Noninvasive Physiology Measurement

Record great wireless data in the lab and in the real world!

BioNomadix Transmitters offer either two of the same signal or a combination of signals. Use with leads, electrodes, transducers, and accessories—including the BioShirt sensing shirt. Wirelessly transmit to a paired Receiver module, Smart Center, or Logger.



- | | |
|------------------|---------------|
| Dual ECG | Pulse+EDA |
| Dual EEG | Resp+ECG |
| Dual EGG | NICO (dZ/Dt) |
| Dual EMG | Accelerometry |
| Dual EOG | Clench-EMG |
| Dual Respiration | Goniometry |
| Dual Skin Temp | Heel/Toe |

Uncompromised participant comfort and freedom of movement for an array of life science applications—Psychophysiology, Exercise Physiology, Biomechanics, Consumer Neuroscience, Cardiology—HRV, Evoked Response, and more!



Record Great Physiology Data Where, When, and How You Want

BioNomadix Receiver Sets

Matched Transmitter/Receiver sets deliver data in a wireless system

- Wireless, wearable physiology anywhere
- Small & unobtrusive
- Reduces cables for greater flexibility
- Record from multiple subjects without interference
- Use with MP160 System or third-party DAQ. Each MP160 records up to 16 channels; also works with multiple MP160 units



BioNomadix Smart Center



ACQKNOWLEDGE[®] for Smart Center



- **Power with Portability**
Groundbreaking portability, in-the-lab or on-the-go. The Smart Center is engineered for powerful data collection, visualization and analysis.
- **Sophisticated and Streamlined**
Set up and start recording in minutes with the easy-to-use auto-configuration wizard, auto-detection of transmitters and analysis recommendations.
- **Compact, Convenient, and Complete**
Flexible system allows subjects to be mobile, remote, and comfortable.

BioNomadix Logger

Wirelessly record great data in the real world

- Get great physiological data while participants live their lives
- Built-in color display, speaker, vibration, voice journal, event markers, alarms, and accelerometer for activity info
- Perfect tool for applications that demand greater degrees of subject freedom and complex experimental design
- Add GPS tracking and synchronize with subject's audio notes



Physiology where, when, and how you want

The BioNomadix® system of wearable wireless devices delivers the freedom to discover the data the researcher desires, in the environment and at the scale of the researcher's choosing, with the quality scientific research demands, and an unparalleled ease of use for both researcher and subjects.



Record great wireless data in the lab and in the real world!

The latest generation BioNomadix Transmitters can operate with the BioNomadix Smart Center, Logger, BioShirt, or matched Receiver to noninvasively record full signal bandwidth physiology data (existing Transmitters require a firmware update for Logger capability).

Dual-Signal Transmitters

- ECG** **Respiration**
- EEG** **Temperature**
- EMG** **Cardiac Output**
- EOG** **Heel & Toe Strike**
- EGG** **Clench Force**
- EDA** **Accelerometer**
- Pulse** **Goniometry**



When transmitted via paired Receiver to an MP160 System, up to 16 channels of BioNomadix® data can be recorded for multi-subject or multi-parameter protocols. The system also works with multiple MP160 systems or third-party data acquisition hardware via an isolated power supply module.

BioNomadix accessory items include transducers, electrode leads, straps, and shirts. The new BN-BIOSHIRT contains a respiration sensor and fabric electrodes to simultaneously acquire Respiration and ECG data from freely roaming participants, and for multi-sensor protocols, the pocketed, stretch-mesh BN-SHIRT comfortably holds multiple devices.



BioNomadix is the perfect tool for applications that demand greater degrees of subject freedom and complex experimental design!


BioNomadix with AcqKnowledge deliver the quality data your research demands

BioNomadix high-fidelity wireless recording and AcqKnowledge software provide a powerful, complete, wireless solution that supports advanced analysis for applications and measurements for a variety of physiological parameters, including: EDA, Heart rate, respiration rate, Heart rate variability (HRV), Respiratory Sinus Arrhythmia (RSA), etc. Combine BioNomadix data with other signals for a comprehensive analysis of the participant's experience—GPS, Eye Tracking, Video, etc.



AcqKnowledge software displays, controls, analyzes, replays, and exports BioNomadix data in one convenient program. Plus, AcqKnowledge provides the power of sophisticated automation and scoring routines for each signal type, and customization options.

New AcqKnowledge for Smart Center Wizard streamlines setup and provides signal optimized analysis options. Start collecting data within minutes and easily store/retrieve multiple protocols for repeatability.



For greater freedom, use BioNomadix wireless physiology with a BIOPAC Research System

Combine the sophistication and performance of BIOPAC data acquisition hardware with the power and flexibility of AcqKnowledge software to customize your acquisition and analysis system for life science research.

MP160 data acquisition system




- High resolution — 16 bit
- High speed — up to 400 kHz aggregate
- Variable sample rates (analog & calculation channels)
- 16 analog inputs and 2 independent analog outputs
- Digital I/O lines (receive/send TTL triggers)
- 16 online calculation channels
- Ethernet connectivity — fast and efficient
- Safety

ACQ KNOWLEDGE

The AcqKnowledge software included with each MP System is a highly interactive user-friendly application with intuitive controls to instantly view, measure, analyze, and transform data.

Perform complex data acquisition, triggering and analysis using simple pull-down menus and dialogs — no need to learn a programming language or new protocol.

- **Analysis Features** — signal averaging, sophisticated pulmonary integration routines, filtering, FFT, histogram, automatic data reduction, template analysis, peak detection features, find rate settings, and an equation generator
- **Remote Monitor** — view subject data on another machine – bedside monitor display. Simplified user interface tracks the welfare of the subject with alarms to warn when signals fall out of range. The system will work on any device that has access to the same IP based network as the MP160 or MP150.



- **Multi-media Support Tools** — videos for analysis, automation routines, hardware setup, and scripting; real-time, searchable user guides and application notes (PDF), as well as extensive online support, knowledge base, and training options. Plus, Module Setup Wizard and QuickStart template files are included to make it even easier to start your experiment.


- **Acquisition Features** — variable sample rates, pause mode, and stimulation design and control. Online analysis settings, filters and transformations provide real-time recording feedback.
- **Rich Display Features** — multiple display modes, advanced grid system, journal facility for note taking, textual event markers, and measurement tools. Mouse-over tool tips (for sample rate, channel rate, measurement results, etc.) help guide application use.
- **Automated Analysis Routines** — save time and standardize interpretation of results with scoring and analysis routines for ECG, HRV, EDA, EMG, EEG, ICG, BP, LVP, Pulmonary Function, and more!
- **BIOPAC Developer** — Customize and automate your analysis routines with new Developer tools including BIOPAC Basic Scripting, Network Data Transfer, or API.

AcqKnowledge software is included with BIOPAC Research Systems and provides comprehensive tools to simplify & standardize advanced analysis.

Standard Analysis Tools

- | | | |
|---------------------|----------------------|-------------------------------------|
| Cycle/Rate Detector | Expression Evaluator | Autoregressive Modeling |
| Digital Filters | Ensemble Averaging | AR Time-Freq Analysis |
| FFT & PSD | Template Analysis | Independent Component Analysis |
| Histogram | Wavelet Analysis | Principal Component Analysis |
| Stim-Response | Nonlinear Modeling | Plotting options - 3D, overlap, X/Y |

Automated Advanced Analysis



- ECG** - Electrocardiography
 - ECG Interval Extraction
 - Heart Rate Variability
 - Chaos Analysis
 - Detect & Classify Heartbeats
- EEG** - Electroencephalography
 - Compute Approximate Entropy
 - Delta Power Analysis
 - Derive Alpha-RMS
 - Derive EEG Frequency Bands
 - EEG Frequency Analysis
 - Remove Artifact - EOG/eye blink
- EMG** - Electromyography
 - Derive Average Rectified EMG
 - Derive Integrated EMG
 - Derive Root Mean Square EMG
 - EMG Frequency & Power Analysis
 - Locate Muscle Activation
- EDA** - Electrodermal Activity
 - Derive Phasic EDA from Tonic
 - Event-related EDA Analysis
 - Locate SCRs
- NICO** - Impedance Cardiography
 - Derive dZ/dt from Raw Z
 - Classify dZ/dt : B, C, X, Y, and O
 - Remove dZ/dt Motion Artifacts
 - ICG Analysis
 - VEPT (electrically participating tissue)
 - PEP Pre-ejection Period
 - Body Surface Area & Ideal Body Weight
- RSP** - Respiration
 - Respiratory Sinus Arrhythmia
- EGG** - Electrogastrography
 - Gastric Wave Analysis
 - Gastric Wave Coupling



TUTORIALS

Watch dozens of tutorial screencasts and learn more online!



Researchers can increase the complexity of experimental design and create complex, real-world scenarios. BioNomadix allows a natural, unhindered environment, which significantly improves the quality of the data and makes it much easier for subjects to achieve peak performance.

The untethering of the subject and elimination of major cables greatly enhances the recording experience and helps to provide a relaxed environment for sensitive study populations that don't like to be tethered.



BioNomadix works extremely well for protocols that demand movement, such as exercise physiology and sports medicine

- BioNomadix provides excellent signal quality with digital transmission and short leads placed close to the signal source.
- The small and unobtrusive nature of the BioNomadix helps to relax the subject and reduce stress and anxiety.
- For multi-subject applications, BioNomadix allows multiple devices to work seamlessly inside an enclosed area.
- The lightweight, comfortable BioShirt transmits ECG and Respiration data and delivers quick setup, great data, and long-term comfort.
- Subjects can freely move around without cables anchoring them to recording devices and the fear of tripping over other subject's cables.
- BioNomadix is small and easy to transport for offsite and home recording.
- The BioNomadix system scales extremely well, from limited-channel applications to multi-sensor studies.

No matter what your application, BioNomadix help provide a more natural, relaxed environment for in lab, real-world, or virtual reality immersion studies



Full-bandwidth Wireless Data

BioNomadix wearable, wireless solutions for life science data record and analyze a variety of physiological parameters in data logging or telemetry modes. Devices incorporate internal, non-distorting, highpass and lowpass filters to provide for high quality amplification of the complete waveform resulting in exceptional quality data. The signal is transmitted via an ultra-low power, 2.4 GHz bi-directional digital RF transmitter. The recording transmitter is extremely safe to use because there is no physical connection to supply mains. Interference with other recording modules or systems is greatly minimized because the module transmitter is completely battery operated and operates totally independently of any other recording grounds or power.

Transmitters

BioNomadix Transmitters offer either two of the same signal or a combination of signals. Use with leads, electrodes, transducers, and accessories. Wirelessly transmit to a paired Receiver module, Smart Center, or Logger.

Dual ECG	BN-ECG2-T	Pulse+EDA	BN-PPGED-T
Dual EEG	BN-EEG2-T	Resp+ECG	BN-RSPEC-T
Dual EGG	BN-EGG2-T	dZ/Dt	BN-NICO-T
Dual EMG	BN-EMG2-T	Accelerometry	BN-ACCL3-T
Dual EOG	BN-EOG2-T	Clench-EMG	BN-DYNEMG-T
Dual Respiration	BN-RSP2-T	Goniometry	BN-GONIO-T
Dual Skin Temp	BN-SKT2-T	Heel/Toe	BN-STRIKE-T

Receiver Sets

BioNomadix Receiver Sets consist of a Transmitter that the subject wears and a paired Receiver module that connects to an MP System—each MP System can record up to 16 channels of data for multi-subject or multi-parameter protocols; also works with multiple MP Systems. For third-party DAQ hardware, add isolated power (IPS100C).

BN-ECG2	BN-EEG2	BN-EGG2	BN-EMG2	BN-EOG2
BN-RSP2	BN-SKT2	BN-PPGED	BN-RSPEC	BN-NICO
BN-ACCL3	BN-DYNEMG	BN-GONIO	BN-STRIKE	

Smart Center

BioNomadix Smart Center is a stand-alone wireless system that includes Transmitters, AcqKnowledge for Smart Center software with streamlined setup and analysis tools, and a carrying case for storage or portability; the Enhanced System adds a Logger and Scripting.

Smart Center Essentials with 2 Transmitters **BN-SMART-ESS2**

3 Transmitters **BN-SMART-ESS3**

Smart Center Enhanced with 3 Transmitters, Logger, Scripting **BN-SMART-ENH**

Logger

The **BioNomadix Logger** records up to 24 hours from up to 3 Transmitters plus its built-in Accelerometer, for up to 10 channels of data. **BN-LOGGER-1/2/3** Stand-alone system with Logger, AcqKnowledge and Transmitter(s); **BN-LOGGER** Logger with AcqKnowledge; **BN-LOGGER-ADD** Logger unit only (for an existing Smart Center or Receiver Set system).

Use with BioNomadix Wireless Transmitters

Dynamometry

BN-CLENCH-XDCR Clench bulb dynamometer measures proportionality of clench force to pressure on the bulb.



EEG Cap System

BN-EEGCAP-SYS Cap system (19 electrodes, 10/20 montage) with accessories. Requires **BN-ADAPT-TP2/TP3**. **BN-CAP-SIZE** Small 50-54 cm, Medium 54-58 cm, Large 58-62 cm



Goniometry*

BN-GON-110-XDCR, BN-GON-150-XDCR Twin-axis transducer with two separate outputs — measure flexion/extension and/or radial/ulnar deviation. Choose 110 mm or 150 mm.
BN-GON-F-XDCR Finger goniometer, 35 mm, measures polycentric joint movement.



Pulse

BN-PULSE-XDCR, BN-PULSEEAR-XDR Infrared emitter and photo-diode transmit changes in infrared reflectance from varying blood flow; finger or ear.



Respiration

BN-RESP-XDCR Measure changes in thoracic or abdominal circumference as subject breathes.



Strike

BN-STRIKE-XDCR Record heel/toe strike activity as subject walks.



Temperature

BN-TEMP-A-XDCR Skin temp - reponse time 1.1 sec. Stainless steel banjo design; tape to body.
BN-TEMP-B-XDCR Fast response 0.6 sec (in air)



Torsiometry*

BN-TOR-110-XDCR, BN-TOR-150-XDCR Torsiometers measure axial rotation about a single plane. Choose 110 mm or 150 mm.



*Adapters

BN-GON/BN-TOR include required adapter(s) to interface with BioNomadix transmitter. To use existing BIOPAC or 3rd-party goniometers/torsiometers, order one **BN-ADAPT-GONIO** per channel.

Leads and Electrodes

Mini-pinch leads two-lead or three-lead in 15 cm, 30 cm, or 45 cm

BN-EL15-LEAD2 **BN-EL30-LEAD2** **BN-EL45-LEAD2**
BN-EL15-LEAD3 **BN-EL30-LEAD3** **BN-EL45-LEAD3**

EDA leads 2 x 15 cm or 2 x 25 cm, use with disposable EDA electrodes

BN-EDA-LEAD2 **BN-EDA25-LEAD2**

NICO leads 50 cm clip leads - 2 x alligator or 4 x mini-pinch

BN-EL50-LEAD2 **BN-EL50-LEAD4**

Touchproof adapters 15 cm, sockets accept TP leads

BN-ADAPT-TP2 (+, -) **BN-ADAPT-TP3** (+, -, GND)

Disposable electrodes vinyl or cloth, pre-gelled or dry **EL500 series**



Shirts

Sensing Shirts **BN-BIOSHIRT-SIZE**



BioShirts provide simultaneous ECG and Respiratory signals from sensors that are fully integrated into the shirt's textile structure; use with BN-RSPEC-T Transmitter and paired Receiver or Logger (separate purchase)—shirt pocket holds transmitter. Comfortable shirts can be worn under clothing. Ideal for exercise and long-term studies. Designed for snug fit.



FXS-FXL Female extra small through extra large
MS-MXXL Male small through extra extra large

Shirt

BN-SHIRT-size (XS, S, M, L, XL availability varies by size)



When multiple transmitters are used, this custom stretch-mesh shirt provides a greater degree of comfort and mounting flexibility. Pockets, zippers and bands help to place and hold transmitters. Select size for compression fit to hold the BioNomadix transmitter and sensors in place; wear as is or under clothing.

Accessories

Straps

BN-STRAP-size (20, 33, 76, 137 cm)

Each system includes a strap for common signal type applications. Additional straps are available to fit wrist, head, leg, chest, etc. (sized in cm: BN-STRAP-20 is 20 cm).

Chargers

BN-BAT-CHRG for Transmitters **BN-LOG-CHRG** for Loggers

Transmitter charger typically provides full charge (70-90 hrs) in one hour; 500 cycles. Logger charger typically reaches full charge (24 hr acquisition, 30 day standby) in 12 hrs.

I/O Interface **BN-SMART-IOCBL**

Access 8 Digital TTL Inputs on the Smart Center; 1 m cable from I/O port to Dsub25 male.

Table 1: BioNomadix Smart Center

Unit Dimensions:	92 mm (L) x 60 mm (W) x 27 mm (D)
Maximum Sample Rate:	2 kHz per channel
Bit Rate:	12 bits per sample
Frequency:	2.4 GHz
Ports:	USB (1), I/O (1)
Antenna:	Yes, removable
Power Source:	Computer USB (cable included)
Range:	10 meters line-of-sight
Transmission:	Wireless—with BioNomadix Transmitters
Maximum Transmitters:	3 Dual-Signal or Tri-Signal Transmitters
Maximum Data Channels:	9, plus derived signals selectable in software
FCC ID IC VCCI:	ZWIBN XR1 9901A-BN XR1 211-128161
Software:	AcqKnowledge 5 for Smart Center
Supported OS:	Windows 10/8.x/7, Mac OS X 10.10-10.12
Carrying Case Dimensions:	34 cm (L) x 28.47 cm (W) x 8.24 cm (D)

Table 2: BioNomadix Logger & Logger Charger

BN-LOGGER	
Transmitter:	Ultra-low power 2.4 GHz bi-directional digital RF transmitter
Built-in Accelerometer:	X, Y, Z- axes; rate 100-400 Hz; Range: 2-16 G
Rate	2 kHz, maximum
RF Reception range:	1 meter (line of sight, approx.)
Memory:	8 GB
Operating time:	24 hours (recording)
Battery:	1800 mAh Lithium-ion
Screen:	Color, 6 cm diagonal
Dimensions:	9.42 cm x 5.76 cm x 2.3 cm
Weight:	121.2 grams
Compliance:	FC, CE, IC, VCCI -FCC Part 15 B FCC ID: ZWIBNXT1, IC: 9901A-BNXT1
BN-LOG-CHARGER	
Charger style:	Integrated USB charger with AC wall adapter
Operating time:	Acquisition 24 hours; stand-by ~30 days
Time to full charge:	~12 hours

Table 3: Common BioNomadix Specs—Receiver Modules & Transmitters

Operational Range:	10 meters (line-of-sight) typical in standard laboratory setups See also: <i>Operational Range and Characteristics</i>
Delay:	Large fixed component (15.6 ms) and small variable component (± 0.5 ms rms)
Temp & Humidity:	Operating Temperature: 5-45° C Humidity: 95% non-condensing
Size & Weight:	Transmitter (approx.): 6 cm x 4 cm x 2 cm 54 grams Receiver (approx.): 4 cm x 11 cm x 19 cm 380 grams
Transmitter:	Type: Ultra-low power, 2.4 GHz bi-directional digital RF transmitter Rate: 2,000 Hz (between transmitter and receiver)
Receiver Power:	Use with an MP Research System or with isolated power supply IPS100C for 3rd-party data acquisition system.
Battery:	BioNomadix transmitters use an L-ion battery: full charge takes ~1 hour to provide maximum operating time.
Transmitter Charger:	A battery charger is included with each module pair. See BN-CHARGER for charge time and recharge cycle details.
Compliance:	FCC, CE, IC, - FCC Part 15 B - FCC ID: receiver: ZWIBN XR1, transmitter: ZWIBN XT1 IC: receiver: 9901A-BN XR1, transmitter: 9901A-BN XT1

Table 4: BioNomadix Dual Biopotential Pairs

BioNomadix	BN-ECG2	BN-EEG2	BN-EGG2	BN-EMG2	BN-EOG2
Signal type:	Dual Channel ECG	Dual Channel EEG	Dual Channel EGG	Dual Channel EMG	Dual Channel EOG
Bandlimits Max:	0.05 Hz to 150 Hz	0.1 Hz to 100 Hz	0.005 Hz to 1.0 Hz	5 Hz to 500 Hz	0.005 Hz to 100 Hz
Factory preset:	1 Hz to 35 Hz	0.5 Hz to 35 Hz	0.005 Hz to 1.0 Hz	10 Hz to 500 Hz	0.005 Hz to 35 Hz
Filter options:	0.05 or 1 Hz HP, 35 or 150 Hz LP	0.1 or 0.5 Hz HP, 35 or 100 Hz LP	0.005 Hz HP, 1 Hz LP	5 or 10 Hz HP, 250 or 500 Hz LP	0.005 or 1 Hz HP, 35 or 100 Hz LP
Alternative signal:	Heart Rate Mode	Delta, Theta, Alpha, Beta		Envelope Detection Mode	Derivative Mode
Notch filter:	50/60 Hz user-controlled switch; typically not required—factory preset OFF. See Appendix for more hardware-specific output options.				
Noise Voltage (shorted inputs):	0.9 μ V rms (bandwidth of 0.05 Hz to 150 Hz)	0.2 μ V rms (bandwidth of 0.10 Hz to 100 Hz)	0.5 μ V rms (bandwidth of 0.005 Hz to 1 Hz)	1.5 μ V rms (bandwidth of 1.0 Hz to 500 Hz)	0.9 μ V rms (bandwidth of 0.005 Hz to 100 Hz)
Input Voltage:	up to 10 mV P-P	up to 2 mV P-P	up to 10 mV P-P	up to 10 mV P-P	up to 10 mV P-P
Output Voltage:	± 10 V (receiver output)				
CMRR	110 dB typical at 50/60Hz; 90dB minimum for ECG, EEG, EMG, and EOG, 100 db minimum for EGG				
CMII & Imped.	CMII: 1000 M Ω (50/60 Hz) Differential Input Impedance: 2 M Ω				
Fixed Gain:	2,000	10,000	2,000	2,000	2,000
Operating Time:	72-90 hours				
Included strap:	137 cm BN-STRAP137	76 cm BN-STRAP76	137 cm BN-STRAP137	33 cm BN-STRAP33	76 cm BN-STRAP76
Size (approx.) & Weight:	Transmitter: 6 cm x 4 cm x 2 cm 54 grams			Receiver: 4 cm x 11 cm x 19 cm 380 grams	
Input:	See BioNomadix electrode lead cable options (BN-ELxx-LEADx). Each biopotential transmitter requires at least one GND. To eliminate redundant biopotential GND, use a 3-lead electrode lead cable for one input (CH A or B) and a 2-lead electrode lead cable for the other input (CH A or B) on each BioNomadix transmitter. Use BN-ADAPT-TP2/3 for Touchproof connections, including BN-EEGCAP-SYS.				

Specifications subject to change without notice. (6/2019)

Table 5: BioNomadix Dual Transducer Pairs

BioNomadix	BN-SKT2	BN-RSP2	BN-GONIO	BN-STRIKE
Signal type:	Dual Channel SKT temp	Dual Channel RSP respiration	Dual Channel Goniometry	Dual Channel Strike Data
Bandlimits Max:	DC to 10 Hz	DC to 10 Hz	DC to 100 Hz	DC to 100 Hz
Factory preset:	DC to 1 Hz	DC to 1 Hz	DC to 10 Hz	DC to 10 Hz
Filter Options:	DC, 0.5 Hz HP, 1 or 10 Hz LP	DC, 0.5 Hz HP, 1 or 10 Hz LP	DC, 3 Hz, 10 Hz, or 100 Hz LP	DC, 3 Hz, 10 Hz, or 100 Hz LP
Notch filter:	50/60 Hz user-controlled switch—factory preset OFF; typically not required. See Appendix for additional hardware-specific output options.			
Resolution:	0.01° C (rms)	FSR/4096; (4.88 mV)	0.01° rotation (rms)	N/A
Signal range:	13 to 51° C	± 10 V (at output)	± 180°	± 10 V (at output)
Output Voltage:	± 10 V (receiver output)			
Operating time:	72-90 hours			
Included strap:	137 cm BN-STRAP-137	137 cm BN-STRAP-137	76 cm BN-STRAP-76 & BN-STRAP-33	33 cm BN-STRAP-33
Input:	BN-TEMP-A/B-XDCR	BN-RESP-XDCR	BN-GON-110-XDCR BN-GON-150-XDCR BN-GON-F-XDCR BN-TOR-100-XDCR BN-TOR-150-XDCR	BN-STRIKE-XDCR

Table 6: BioNomadix Accelerometer

BioNomadix	BN-ACCL3
Signal type:	G (X, Y, Z)
Signal range:	Selectable: ±2, ±4, ±8 or ±16 G
Bandlimits Max:	±2, ±4, ±8 or ±16 G
Factory preset:	± 16 G at 400 Hz LP
Filter Options:	DC to 3.13 Hz LP up to 400 Hz LP (in power of 2 steps)
Alternative:	Tap Event Mark Mode (replaces G)
Resolution:	X: 5 mg rms, Y: 6 mg rms, Z: 9 mg (rms) (±2 G scale at 400 Hz LP)
Output Voltage:	±10 V (receiver output)
Operating time:	72-90 hours
Included strap:	33 cm - BN-STRAP33
Input:	Attach BioNomadix transmitter to subject—no additional hardware input required; sensor is internal to transmitter

Table 7: BioNomadix Signal Combo Pairs

BioNomadix	BN-RSPEC	BN-PPGED	BN-NICO	BN-DYNEMG
Signal type:	RSP plus ECG	PPG plus EDA	Z and dZ/dt	Dynamometry plus EMG
Bandlimits Max:	Respiration (CH A): see BN-RSP2 spec	DC to 10 Hz; 0.5 Hz to 3 Hz	Both: DC to 50 Hz	DYN: DC 100 Hz
Factory preset:	ECG (CH B): see BN-ECG2 spec	DC to 3 Hz	Both: DC to 50 Hz	DYN: DC to 10 Hz
Filter Options:		DC, 0.5 Hz HP, 3 or 10 Hz LP	DC, 1, 3, 5, 50 Hz LP	DYN: DC, 3 Hz, 10 Hz, or 100 Hz LP
		1 Hz LP	Excitation: Type: Alternating current sine wave, Current: 1 mA, rms Frequency: 50 kHz	EMG: see BN-EMG2 specs
Notch filter:	50/60 Hz user-controlled switch; typically not required—factory preset OFF. See Appendix for additional hardware-specific output options.			
Resolution:	see BN-RSP2 and BN-ECG2 specs	PPG: FSR/4096; (4.88 mV) EDA: 0.012 µS (min step)	Z: nominally -0.05 Ω (rms) at 10 Hz BW dZ/dt: -0.0075 Ω/sec (rms) at 10 Hz BW	DYN: 35 micro kg-f/cm2 (0.0005 psi) (rms) EMG: see BN-EMG specs
Signal range:	see BN-RSP2 and BN-ECG2 specs	PPG: ±10 V (at output) EDA: 0 to 50 µS; excitation: 0.5 V constant V	Z: 5 to 100 Ω (mag) dZ/dt: ±10 Ω/sec	DYN: 0-1.055 kg-f/cm2 EMG: up to 10 mV P-P
Output Voltage:	± 10 V (receiver output)			
Operating time:	72-90 hours	24 hours	24 hours	75 hours
Included strap:	137 cm BN-STRAP137	33 cm BN-STRAP33	137 cm BN-STRAP137	33 cm BN-STRAP-33
Input:	CH A: BN-RESP-XDCR CH B: BN-ELxx-LEAD3	CH A: BN-PULSE-XDCR or BN-PULSEEAR-XDR CH B: BN-EDA LEAD2 or BN-EDA25-LEAD2	2 x BN-EL50-LEAD4 or 2 x BN-EL50-LEAD2	CH A: BN-CLENCH-XDCR CH B: BN-ELxx-LEAD3



Wireless, Wearable Physiology
Life Science Monitoring Devices
for Live or Logged Data



Two-Channel Transmitters

- Dual EC
- Dual EEG
- Dual EOG
- Dual EMG
- Dual EGG
- Dual Temperature
- Dual Respiration
- Respiration with ECG
- Electrodermal Activity with Pulse
- Impedance Cardiography
- Heel/Toe Strike
- Dynamometry with EMG
- Goniometry
- Torsiometry
- Accelerometry

www.biopac.com info@biopac.com



42 Aero Camino
Goleta, CA 93117