



BIOPAC
Systems, Inc.

Physiology Lessons
for use with the
Biopac Student Lab

PC under Windows® 98SE, Me, 2000 Pro
or Macintosh® OS 8.6-9.1

Manual Revision
PL3.6.7-ML3.0.7/061903

Richard Pflanzler, Ph.D.
Associate Professor
Indiana University School of Medicine
Purdue University School of Science

J.C. Uyehara, Ph.D.
Biologist
BIOPAC Systems, Inc.

William McMullen
Vice President
BIOPAC Systems, Inc.

BIOPAC Systems, Inc.




42 Aero Camino, Goleta, CA 93117
(805) 685-0066, Fax (805) 685-0067

Email: info@biopac.com

Web Site: <http://www.biopac.com>

Lesson 5 Data Report

ELECTROCARDIOGRAPHY I
Components of the ECG



Lesson 5

ELECTROCARDIOGRAPHY I

ECG I

DATA REPORT

Student's Name: _____

Lab Section: _____

Date: _____

I. Data and Calculations

Subject Profile

Name _____

Height _____

Age _____

Weight _____

Gender: Male / Female

A. Supine, Resting, Regular Breathing (using *Segment 1* data)

Complete the following tables with the lesson data indicated, and calculate the Mean and Range as appropriate.

Table 5.3

Measurement	From Channel	Cardiac Cycle			Mean	Range
		1	2	3		
ΔT	CH 2					
BPM	CH 2					

Table 5.4

ECG Component	Duration ΔT [CH 2]				Amplitude (mV) Δ [CH 2]			
	Cycle 1	Cycle 2	Cycle 3	Mean	Cycle 1	Cycle 2	Cycle 3	Mean
	P wave							
PR interval								
PR segment								
QRS complex								
QT interval								
ST segment								
T wave								

Table 5.5

Ventricular Readings	CH 2 ΔT			
	Cycle 1	Cycle 2	Cycle 3	Mean
QT Interval <i>(corresponds to Ventricular Systole)</i>				
End of T wave to subsequent R wave <i>(corresponds to Ventricular Diastole)</i>				

B. Seated, deep breathing

Table 5.6

Rhythm	CH. #	Cycle 1	Cycle 2	Cycle 3	Mean
<i>Inspiration</i>					
ΔT	CH 2				
BPM	CH 2				
<i>Expiration</i>					
ΔT	CH 2				
BPM	CH 2				

C. Sitting

Table 5.7

Heart Rate	CH. #	Cycle 1	Cycle 2	Cycle 3	Mean
ΔT	CH 2				
BPM	CH 2				

D. After Exercise

Table 5.8

Ventricular Readings	CH 2 ΔT			
	Cycle 1	Cycle 2	Cycle 3	Mean
QT Interval <i>(corresponds to Ventricular Systole)</i>				
End of T wave to subsequent R wave <i>(corresponds to Ventricular Diastole)</i>				

Measurement	Mean	Range
<i>Supine, regular breathing</i>		
Ventricular systole	_____	_____
Ventricular diastole	_____	_____
<i>After Exercise</i>		
Ventricular systole	_____	_____
Ventricular diastole	_____	_____

What changes occurred in the duration of systole and diastole between resting and post-exercise?

G. Review your Data

1. Is there always one P wave for every QRS complex? Yes No
2. Describe the P and T wave shapes: _____
3. Do the wave durations and amplitudes for all subjects fall within the normal ranges listed in Table 5.2? Yes No
4. Do the ST-segments mainly measure between -0.1 mV and 0.1 mV? Yes No
5. Is there baseline “drift” in the recording? Yes No
6. Is there baseline “noise” in the recording? Yes No

End of Lesson 5 Data Report